





**SUPER  
LUX**  
MOULD CHART

## Upper Anteriors



50	7,0	37,7	8,5
1/44	7,6	42,0	9,5
1/13	7,6	43,4	9,5
1/47	8,2	46,3	9,8
1/17	8,6	47,3	10,3
1/40	8,7	48,0	11,0
1/49	8,7	48,5	11,0





## Upper Anteriors



1/30	6,9	38,4	8,0
56	7,6	42,9	8,1
58	7,8	44,0	9,0
52	8,0	44,5	10,0
1/48	8,5	46,2	10,0
1/37	8,5	46,5	10,1
59	8,6	47,0	10,0

### STRUCTURE

GESICHSFORM  
 STRUCTURE  
 ESTRUCTURAS  
 CONFORMAZIONE  
 UDFORMING  
 MALLI  
 ΜΟΡΦΗ  
 CONFORMATIE  
 CONFORMAÇÃO  
 FORM



## Upper Anteriors



## Lower Anteriors




AVAILABLE SHADES

**A1 - D4**  
**1A - 6D**

## MESIO-DISTAL PROFILE

PROFIL MESIO-DISTAL  
 PROFIL MESIO-DISTAL  
 PERFILE MESIO-DISTAL  
 PROFILO MESIO-DISTALE  
 MESIO-DISTALPROFIL  
 MESIO-DISTAALINENPROFILI  
 ΟΨΗ MESIO-DISTAL  
 PROFIEL MESIO-DISTAL  
 PERFIL MESIO-DISTAL  
 MESIO-DISTALPROFIL




CURVED

FLAT

## LATERAL PROFILE



PROFIL LATERAL  
 PROFIL LATERAL  
 PERFILE LATERAL  
 PROFILO LATERALE  
 SIDEPROFIL  
 SIVUPROFILI  
 Πλάγια όψη  
 ZIJPROFIEL  
 PERFIL LATERAL  
 SIDOPROFIL



CURVED

FLAT

## L-cusp Posteriors



## SUPER LUX

- T3011 Anteriors Set A1-D4 x 6
- T3021 Anteriors Set 1A-6D x 6
- T3012 Posteriors Set A1-D4 x 8
- T3022 Posteriors Set 1A-6D x 8
- T3019 Complete Set A1-D4 x 28
- T3029 Complete Set 1A-6D x 28



## CHOICE SYSTEM

in relation to the size

AUSWAHLSYSTEM NACH DEM KÖRPERBAU  
 SYSTEME DE CHOIX EN RELATION A LA TAILLE  
 METODO DE ELECCION EN BASE A LA COMPLEXION FISICA  
 MISURAZIONE IN BASE ALLA CORPORATURA  
 MÅLING EFTER KROPPSBYGNINGEN  
 MITOITUS RUUMINRAKENTEEN PERUSTEELLA  
 ΜΕΤΡΗΣΗ ΒΑΣΕΙ ΣΩΜΑΤΙΚΗΣ ΔΙΑΠΛΑΣΗΣ  
 METINGEN OP BASIS VAN LICHAAMSSPECIFIEKE VORMEN  
 MEDIÇÃO SEGUNDO A ESTATURA  
 MÅTNING PÅ BASIS AV KROPPSBYGGNAD




L

Leptosonal



P

Pyknic



A

Athletic

## COMBINATION CHART

ARTIKULATIONS-TABELLE  
 TABLEAU D'ARTICULATION  
 TABLA DE ARTICULACION  
 TABELLA DI ARTICOLAZIONE  
 OVERSIGHT OVER ARTIKULATION  
 PARENTATAULUKKO  
 ΠΙΝΑΚΑΣ ΑΡΘΡΩΣΕΩΝ  
 ARTICULATIETABEL  
 TABELA DE ARTICULAÇÃO  
 TABELL ÖVER LEDRÖRELSE

Anteriors		Posteriors	
Upper	Lower	L-cusp	N-cusp
50	0/3	1/60	
1/44	0/4 - 0/5	1/72	76N
1/13	0/5	1/74	76N
1/47	0/8 - 0/6	1/72	77N
1/17	0/8 - 0/10		77N
1/40	0/10 - 0/11		79N
1/49	0/10 - 0/11		79N
1/30	0/3	1/60	70N
56	0/0 - 0/4	1/65	76N
58	0/6 - 0/53	1/72	76N
52	0/53 - 0/6		77N
1/48	0/6	1/72	77N
1/37	0/8 - 0/10		79N
59	0/53 - 0/8	1/72	77N
1/20	0/3	1/60	70N
1/32	0/0 - 0/4	1/65	70N
1/35	0/6 - 0/53	1/72	77N
1/22	0/5 - 0/53	1/74	76N
1/25	0/8 - 0/53		77N
53	0/53 - 0/6	1/72	76N
1/27	0/10 - 0/11		79N
62	0/11		79N

## N-cusp Posteriors



## SUPER LUX SHADE GUIDE



